

## Xela Rederm

- 1 Achieve glowing skin and reduction in superficial wrinkles
- 2 Improves skin pigmentation
- 3 Improved skin tone and reduce dullness of skin



### IMPROVED HYALUAL PURIFICATION TECHNOLOGY

- **Double Synergic Effect Technology**  
The maximum synergistic result of hyaluronic and succinic acid)
- **Monodispersity of HA molecules** <2.5 million Da
- **Enhanced Buffer Protection**

#### TREATMENT CYCLE


3 - 4 sessions at an interval of 2 weeks.




### HOW TO VERIFY PRODUCTS



1. Spot our JLC Global QR Sticker on the box



2. Before treatment, scan the QR code with your phone camera



3. Authenticity verification


exclusively distributed by

**jlc.global**

science driven solutions

🌐 [www.jlcglobal.com.my](http://www.jlcglobal.com.my)


📱 @jlc.global



**XELA  
REDERM**

**SKIN REDERMALIZATION  
THE BEST SOLUTION  
FOR IMPROVING SKIN PIGMENTATION**

# Redermalization




## Hyaluronic Acid

- Activates anti-inflammatory mediators
- Restores and maintains water balance
- Stimulates the synthesis of basic structural proteins
- Provides mechanical support for tissues

## Succinic Acid

- Normalizes microcirculation & activates lymphatic drainage
- Regulates cellular energy
- exchange
- Contains antioxidant
- properties
- Prevents damage of cell
- membrane
- Stimulates the synthesis of basic structural proteins

### Mechanisms of Action



## Xela Rederm products based on HA Concentrations



### Xela Rederm 1.1%

HA 11 mg/ml  
SA 16 mg/ml



### Xela Rederm 1.8%

HA 18 mg/ml  
SA 16 mg/ml



### Xela Rederm 2.2%

HA 22 mg/ml  
SA 16 mg/ml

## TREATMENT RESULTS



\*Results may vary. Unretouched photos were taken before and after the course of procedures